



Commercial Quality Quick Start Guide

1. Battery Operated Machine
2. Charging Cord
3. Reflexology Footpad
4. 4 Rounded Attachments
5. 4 Pointed Attachments
6. 4 Circulation Rings
7. Remote Control

Unbox your machine and make sure that all parts are included.

If anything is missing Innergizerplus@gmail.com call/text (813) 507-4429

Included in your box should be the following:

Before using the machine plug the charger into the power outlet and into the machine and charge fully. During charging the red light will remain lit, the three green lights will light up sequentially as the machine becomes more charged. Once all three of the green lights are lit, the machine will be fully charged and ready to use. Rotate the dial on the base of the machine to turn on the machine. The farther you rotate the dial, the faster the machine will oscillate and the more intense the massage will be. If the machine does not start, the first step should be to make sure the machine is fully charged and unplugged from the charger.

To activate the wireless remote control make sure the machine is fully charged.

Turn the dial until it clicks then the lights will flash sequentially and the remote will be ready to pair by simply pushing the red button. Simply use the +/- buttons to adjust the speed to stop the machine hit the red button.

What makes the Professional InnergizerGo Massager, Superior to others available?

The Professional Massager is 8 Massagers in one.

1. Lifetime Warranty
2. The Innergizer Go is powered by a rechargeable Battery. More Power without a cord.
3. Use the Pad without Attachments (this is the pad attached to the machine, with 4 screw holes); apply it anywhere on the body, even the face.
4. Use the Machine with Rounded Attachments; apply over thicker muscle or skin area to give a general massage.
5. Reflexology Foot Pad is perfect for use on your feet. Simply place over the four rounded attachments on the machine.
6. Use the Machine with the Deep Muscle Attachments; apply on the thickest tissue to work our stress.
7. Use the Machine with the Rounded Circulatory Attachments (caps on deep muscle); apply anywhere on the body for stimulation in the blood circulation.
8. Control the intensity of the machine with either the speed dial on the machine
9. Use the Machine on yourself, simply place it in a chair and lean back, your body weight will hold the machine in place, or on someone else can manipulate the massager on your back for you.



Commercial Quality - Lifetime Warranty

What makes our Professional Innergizer Go Massager, better and different from others on the Market?

The Professional Massager is 8 Massagers in one.

1. LIFETIME WARRANTY
2. The **Innergizer** Go is powered by a rechargeable Battery. More power without a cord.
3. Use the Pad without Attachments (this is the pad attached to the machine, with 4 screw holes); apply it anywhere on the body, even the face.
4. Use the Machine with Rounded Attachments; apply over thicker muscle or skin area to give a general massage.
5. Use the Machine with Reflexology Foot Pad (Rubber foot pad with knobby points); a massage that targets the base of the sciatic nerve and Plantar Fasciitis.
6. Use the Machine with the Deep Muscle Attachments; apply on the thickest tissue to work our stress.
7. Use the Machine with the Rounded Circulatory Attachments (caps on deep muscle); apply anywhere on the body for stimulation in the blood circulation.
8. Use the Machine with the Pyramid Attachments (caps on deep muscle); apply lotions or oils onto the skin, and allow to massage into the skin.
9. Control the intensity of the machine with either the speed dial on the machine or with the optional Phone App (currently only available in Android Store, Apple coming soon). Note; When you turn on the machine with the machine dial, you must turn it off with the machine dial. (Vice Versa for the Phone App) Whichever is turned on first dominates the machine. The other one will not work.
10. Use the Machine on yourself, simply place it in a chair and lean back, your body weight will hold the machine in place, or on someone else

Versatility One Machine can be used anywhere on your body. The intensity of the massage is determined by the pressure and speed you apply.

Power The **Innergizer** Go is an overbuilt Machine, and comes with a Lifetime Warranty. The powerful motor is unlike any other machine on the market. Built to Last with Oscillating Speed from 1300 - 6000 RPM's. Others on the market top out at 3300 RMPs

Attachments Part of our versatility is the assortment of attachments that are included with the Machine.

When you open the box the machine is packaged with the four Rounded Attachments screwed into the top, steel screws with brass fitting, for great quality. This is one of our most popular ways to use the machine. Many people use these four attachments all over the body, specifically on their back, arms and legs. They offer a semi-intense massage that can be varied by the pressure applied to the machine, and/or speed of the machine, and are one of the most popular attachments.

On top of the Rounded Attachment is the Reflexology Foot Pad. This pad is most effectively used on the feet, it is great on the feet for Neuropathy, Plantar Fasciitis and many other conditions. Simply place Rejuvenator Go handle down on the Floor or on a pillow. Place the Reflexology pad over the rounded attachments, or use the Deep Muscle attachments with the Circulation Rings on the machine and place your feet on the massager. While massaging your feet, you will soon discover that the increased circulation and Massage action will be felt all the way up your legs, after using the machine for about 4 -5 minutes you will find that the burst of energy will last for about 45 minutes. Many people who have diminished feeling in their feet and lower legs, have increased sensation after using the machine on a regular basis.

For a more intense Massage you can use the Deep Muscle Attachments in place of the rounded attachments. These attachments are pointed or cone shaped and are able to massage deeper into the muscle for immediate relief. Be careful of applying too much pressure when using these attachments, they are very deep penetrating they should not cause discomfort.

The Circulation Attachments fit as caps on the deep muscle attachments for a different sensation of massage.

How to care for your Innergizer Go Massager:

To store– Place the massager in a safe, dry, and cool area. Avoid contact with sharp edges or pointed objects, anything that may cause damage to the machine. Do not hang or move the unit by its electrical cord.

To clean– Unplug the machine from the power source, remember this whenever handling electronics such as this massager. Allow an hour of cooling after the use of the machine before cleaning. Do not place directly into liquid when cleaning. To clean the surface, run a lightly damp cloth or sponge with detergent or soap, over the surface of the machine. Do not use abrasive cleaners or solvents on this machine.

Important Safety Precautions:

Read all instructions and safety precautions before use.

DANGER-

- To reduce the risk of electric shock, always unplug the appliance from the electrical outlet immediately after charging and before cleaning.
- Never charge the product if it has a damaged cord or plug on the charging cable.
- Do not use while bathing, showering or near any other water source.
- Do not reach for the massager if it has fallen into water, drain water before removing the Rejuvenator Go. Do not use after removing it, for there may be water droplets in the electrical system.
- We recommend using surge protector when charging. Only attach to charge the battery do not leave the **Innergizer Go** plugged in after charging is complete.

WARNING-

- To reduce the risk of burns, fires, electrical shocks, or injury to person. An appliance should never be left unattended when plugged in. Unplug from the electrical outlet when not in use, and before detaching or re-attaching parts.
- Do not operate under blanket or pillow. Excessive heat can occur and fire, electrical shock, or injury may happen.
- Close supervision is necessary when used on, or near children/disabled persons. Do not allow children to operate.
- Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never operate this appliance if it is not working properly, if it has been dropped, damaged or dropped into water. Return the appliance to a service

QUESTIONS?

InnergizerPlus@gmail.com
InnergizerPlus.com

Call/Text 813/507-4429



**Innergizer
Professional Massager**



Deep Muscle



**Circulation
Attachment - Cap**



**Massager on Chair,
held in place by body**



**Upper Back/
Shoulder**



Feet



Forearm



Hand



Calf